Coordinate transformation of correlated images

This example shows how to apply a coordinate transformation previously obtained by correlating a 3D (such as confocal) and a 2D (such as ion beam) image to transform the correlated images. In this way, the confocal image is transformed to the ion beam coordinate system, so that the transformed confocal image "looks like" the ion beam image, or the ion beam image is transformed to the confocal system.

The transformation is read from the output files previously generated by running 3DCT. To make it easier to visualize the transformation, the images used to establish the correlation were artificially generated.

The image transformations shown here are performed using Pyto package.

In [1]:

```
import sys
import pickle
from copy import copy, deepcopy
import numpy as np
import tifffile as tf
import pyto
from pyto.geometry.rigid 3d import Rigid3D
# matplotlib
%matplotlib inline
import matplotlib as mpl
import matplotlib.pyplot as plt
# for convenience
%load ext autoreload
%autoreload 2
%autosave 0
sys.version
```

```
Autosave disabled
Out[1]:
'3.7.3 (default, Oct 7 2019, 12:56:13) \n[GCC 8.3.0]'
```

Small example - Transformation from 3DCT

Read images, markers and transform parameters

```
In [2]:
```

```
# File names
confocal_file_name = 'small_cf.tif'
ion_beam_file_name = 'small_ib.tif'
results_file_name = '2020-04-24_15-03-41_correlation.txt'
```

In [3]:

```
# Rad images and adjust axes
cf_data = tf.imread(confocal_file_name).transpose()
ib_data = tf.imread(ion_beam_file_name)
if len(ib_data.shape) == 3:
    ib_data = ib_data[:,:,0]
ib_data = ib_data.transpose()
```

In [4]:

```
# read marker positions
results = np.loadtxt(results_file_name)
cf_markers = results[:,0:3]
cf_markers_transformed = results[:,3:6]
ib_markers = results[:,6:8]
transform_error = results[:,8:10]
```


In [5]:

```
# Read parameters from 3DCT results
with open(results file name, 'r') as results fd:
    for line in results fd:
        if line.startswith('# - rotation'):
            phi, psi, theta = np.fromstring(line.split(': [')[1].split(']')[0],
sep=',')
            continue
        if line.startswith('# - scale'):
            scale = float(line.split('= ')[1])
            continue
        if line.startswith('# - translation for rotation around [0.0.0]'):
            trans = np.fromstring(line.split('= [')[1].split(']')[0], sep=',')
            continue
        if line.startswith('#
                               translation for rotation around'):
            origin = np.fromstring(line.split('[')[1].split(']')[0], sep=',')
            trans or = np.fromstring(line.split('= [')[1].split(']')[0], sep=','
)
            break
print("Please check that the following parameters are the same as those in the 3
DCT results file:")
print("phi, theta, psi = {}".format((phi, theta, psi)))
print("scale = {}".format(scale))
print("translation around [0, 0, 0] = {}".format(trans))
print("translation around {} = {}".format(origin, trans or))
Please check that the following parameters are the same as those in
```

```
Please check that the following parameters are the same as those in the 3DCT results file: phi, theta, psi = (-143.723, 44.996, 115.579) scale = 0.5 translation around [0, 0, 0] = [-15.361 21.287 0.] translation around [55.55.55.] = [1.75 8.594 -35.247]
```


In [6]:

```
# show max projection confocal
plt.imshow(cf_data.max(axis=2).transpose(), cmap='Greys', origin='lower')
plt.xlabel('x')
plt.ylabel('y')
plt.title('Confocal image (max projection)');
```



In [7]:

```
# Show ion beam image
plt.imshow(ib_data.transpose(), cmap='Greys', origin='lower')
plt.xlabel('x')
plt.ylabel('y')
plt.title('Ion beam image');
```



Transform confocal image to ion beam system


In [8]:

```
# Make transformation from confocal to ion beam (based on 3DCT correlation results)
euler_rad = np.array([phi, theta, psi]) * np.pi / 180
q = Rigid3D.make_r_euler(angles=euler_rad)
tdct = Rigid3D(q=q, scale=scale, d=trans)

# reposition so that transformed markers z > 0
cf_markers_tf = tdct.transform(x=cf_markers, xy_axes='point_dim')
tdct_z = deepcopy(tdct)
tdct_z.d[2] += -np.floor(cf_markers_tf[:,2]).min()
tdct_z.transform(x=cf_markers, xy_axes='point_dim');
```

In [9]:

```
# Transform confocal image
cf_data_tf = tdct_z.transformArray(array=cf_data)
plt.imshow(cf_data_tf.max(axis=2).transpose(), cmap='Greys', origin='lower')
plt.xlabel('x')
plt.ylabel('y')
plt.title("Transformed confocal (max projection)");
```




In [10]:

```
# Transform confocal image with position adjustments
# Usage:
 - Run once with arbitrary shifts and shape
  - Adjust shifts so that all transformed marker coordinates are >0
   - Adjust shape (x, y \text{ and } z) so that x, y \text{ and } z coordinates are smaller
     than the corresponding shape
  - At this point all markers should be visible on the transformed image
#
   - Furher adjust shifts and shape if other parts of the image
     need to be shown on the transformed image
# set shift (reposition) for the transformed image
cf data tf shift = [10, -10, 0]
# set shape of transformed image
cf data tf shape = [60,70,30]
# reposition
tdct local = deepcopy(tdct)
tdct local.d = tdct z.d + cf data tf shift
# transform confocal markers to ion beam system
cf markers tf local = tdct local.transform(x=cf markers, xy axes='point dim')
print("Transformed confocal marker coordinates: \n{}".format(cf markers tf local
))
# transform image
cf data tf = tdct local.transformArray(array=cf_data, shape=cf_data_tf_shape)
plt.imshow(cf data tf.max(axis=2).transpose(), cmap='Greys', origin='lower')
plt.xlabel('x')
plt.ylabel('y')
plt.title("Transformed and shifted confocal (max projection)");
```

```
Transformed confocal marker coordinates: [[45.43334198 5.61596837 0.36918178] [14.61826914 5.58839888 1.39585146] [14.49177391 36.5618814 3.5739021 ] [45.56801854 36.66215395 2.37929325]]
```

Transformed and shifted confocal (max projection)



In [11]:

```
# Save transformed image in tiff format

# transformed image name
cf_tf_file_name = 'small_cf_transformed.tif'

# image data
data = cf_data_tf

# pixel size in um
pixel_um = 0.16

# write tiff
cf_meta = {"PixelSize": str(pixel_um), "FocusStepSize": str(pixel_um)}
cf_res = (10000/pixel_um, 10000/pixel_um, 'CENTIMETER')
tf.imwrite(
    cf_tf_file_name, np.swapaxes(data, 0, 2), metadata=cf_meta, resolution=cf_re
s)
```

In [12]:

```
# Save transformed image in mrc format

# transformed image name
cf_tf_file_name = 'small_cf_transformed.mrc'

# image data
data = cf_data_tf

# pixel size in um
pixel_um = 0.16


# write mrc
cf_tf_image = pyto.core.Image(data=data)
cf_tf_image.write(file=cf_tf_file_name, pixel=1000*pixel_um)
```

Out[12]:

<_io.BufferedWriter name='small_cf_transformed.mrc'>

In [14]:

```
# Optional: Transform confocal image where the transformation is
# defined around the specified center.
# Note: The translation parameters shown here are different form
# those in the corelation results file, because in the latter case
# only rotation center is shifted, but scaling center remains at
# 0 (coordinate system center), while here both rotation and scaling
# centers are shifted.
# specify coordinates of the center
center = (55, 50, 45)
# make transformation and transform markers
tdct or = tdct.shift center(center=center)
cf markers tf or = tdct or.transform(x=cf markers, center=center, xy axes='point
_dim')
# reposition so that transformed markers z > 0
tdct or z = deepcopy(tdct or)
tdct or z.d[2] += -np.floor(cf markers tf or[:,2]).min()
tdct or z.transform(x=cf markers, xy axes='point dim')
print("Parameters of the transformation around the specified center:")
(phi, theta, psi) = tdct_or_z.extract_euler(r=tdct_or_z.q) * 180 / np.pi
print("phi, theta, psi = {}".format((phi, theta, psi)))
print("scale = {}".format(tdct or z.s scalar))
print("translation around {} = {}".format(center, tdct or z.d))
# Transform confocal image and show max projection
cf data tf = tdct or z.transformArray(array=cf data, center=center)
plt.imshow(cf data tf.max(axis=2).transpose(), cmap='Greys', origin='lower')
plt.xlabel('x')
plt.ylabel('y')
plt.title("Transformed confocal (max projection)");
```



Transform ion beam image to confocal system

```
In [15]:
```

```
# Make transformation from ion beam to confocal image
tdct_inverse = tdct.inverse()
```

In [16]:

```
# Convert ion beam marker positions to 3D
ib_markers_3d = np.concatenate((ib_markers, np.zeros((ib_markers.shape[0],1))),
axis=1)
```

In [17]:

```
# Convert ion beam image to 3D

# Set the z-position where the old (2D) ion beam image is
# located within the new (3D) image. This postion should be away from 0
ib_z_position = 5


# convert
ib_data_3d = np.zeros((ib_data.shape[0], ib_data.shape[1], 2*ib_z_position), dty
pe=ib_data.dtype)
ib_data_3d[:,:,ib_z_position] = ib_data
```

In [18]:

```
# Transform ion beam image and show max projection

# set shape of transformed image
ib_data_tf_shape = [90,80,70]

# transform
ib_data_tf = tdct_inverse.transformArray(array=ib_data_3d, shape=ib_data_tf_shape)
plt.imshow(ib_data_tf.max(axis=2).transpose(), cmap='Greys', origin='lower')
plt.xlabel('x')
plt.ylabel('y')
plt.title("Transformed confocal (max projection)");
```




In [19]:

```
# Transform ion beam image with position adjustments
# Usage:
# - Run once with arbitrary shifts and shape
  - Adjust shifts so that all transformed marker coordinates are >0
   - Adjust shape (x, y \text{ and } z) so that x, y \text{ and } z coordinates are smaller
     than the corresponding shape
  - At this point all markers should be visible on the transformed image
#
  - Furher adjust shifts and shape if other parts of the image
     need to be shown on the transformed image
# set shift (reposition) for the transformed image
ib data tf shift = [20, 30, 0]
# set shape of transformed image
ib data tf shape = [120, 140, 100]
# reposition
tdct inverse local = deepcopy(tdct inverse)
tdct inverse local.d = tdct inverse.d + ib data tf shift
# transform ion beam markers markers back to confocal system
ib markers tf local = tdct inverse local.transform(x=ib markers 3d, xy axes='poi
nt dim')
print("Transformed ion beam marker coordinates: \n{}".format(ib_markers_tf_local)
))
# transform image
ib data tf = tdct inverse local.transformArray(array=ib data 3d, shape=ib data t
f shape)
plt.imshow(ib data tf.max(axis=2).transpose(), cmap='Greys', origin='lower')
plt.xlabel('x')
plt.ylabel('y')
plt.title("Transformed confocal (max projection)");
```

```
Transformed ion beam marker coordinates:
[[100.04163704 47.46071089 61.4274567]
[ 55.12025366 31.43912628 21.93740045]
[ 21.28094031 79.7265603 40.8400719]
[ 66.07116356 95.9353055 80.40339432]]
```





```
In [20]:
```

```
# Save transformed ion beam image in tiff format

# transformed image name
cf_tf_file_name = 'small_ib_transformed.tif'

# image data
data = ib_data_tf

# pixel size in um
pixel_um = 0.25

# write tiff
cf_meta = {"PixelSize": str(pixel_um), "FocusStepSize": str(pixel_um)}
cf_res = (10000/pixel_um, 10000/pixel_um, 'CENTIMETER')
tf.imwrite(
    cf_tf_file_name, np.swapaxes(data, 0, 2), metadata=cf_meta, resolution=cf_re
s)
```

In [21]:

```
# Save transformed ion beam image in mrc format

# transformed image name
cf_tf_file_name = 'small_ib_transformed.mrc'

# image data
data = ib_data_tf

# pixel size in um
pixel_um = 0.25

# write mrc
cf_tf_image = pyto.core.Image(data=data)
cf_tf_image.write(file=cf_tf_file_name, pixel=1000*pixel_um, dataType='int16')

Out[21]:
<_io.BufferedWriter name='small_ib_transformed.mrc'>
```

In []: